

## Manual Wheelchair from Europe



**Standard STEEL WHEELCHAIR**

Art. No.VCWK43 **PRICE: 30% below price from the market**

Art . No. VCWK43L **PRICE: 30% below price from the market**

### **PRODUCT CHARACTERISTICS:**

- foldable, steel frame (double cross construction),
- nylon black upholstery,
- soft cushion (VCWK42),
- lifted armrests,
- lifted and removable footrests of adjustable length,
- safety belt,
- calf pad,
- pneumatic quick release rear wheels 24",
- aluminium handrims,
- full front wheels 200x45,
- rear pocket,
- available also without seat cushion (VCWK43),
- frame colours: platinum, blue and green.

### **USE:**

General use wheelchair, at home and for strolls.

	16'	18'
Seat width	42 cm	45 cm
Seat depth	45 cm	45 cm
Wheelchair height	89 cm	89 cm
Wheelchair width	57 cm	63 cm
Wheelchair length	107 cm	107 cm
Width of folded wheelchair	31 cm	31 cm
Weight	18 kg	18 kg
Maximum capacity	100 kg	100 kg



**Standard STEEL WHEELCHAIR**

---

Art. # VCWK42 **PRICE: 30% below price from the market**

---

Art. # VCWK42L **PRICE: 30% below price from the market**

**PRODUCT CHARACTERISTICS:**

- foldable, steel frame (double cross construction),
- nylon black upholstery,
- soft cushion (VCWK42L),
- safety belt,
- lifted armrests,
- rear pocket,
- lifted and removable footrests,
- calf pad,
- full front wheels,
- pneumatic quick release rear wheels 24",
- steel brakes with pressure adjustment,
- varnished steel handrims,
- available also without seat cushion (VCWK43L),
- frame colour: grey.

**USE:**

General use wheelchair, at home and for strolls.

	16'	18'
<b>Seat width</b>	41 cm	45 cm
<b>Seat depth</b>	40 cm	40 cm
<b>Wheelchair height</b>	96 cm	96 cm
<b>Wheelchair width</b>	57 cm	63 cm
<b>Wheelchair length</b>	79 cm	83 cm
<b>Width of folded wheelchair</b>	28 cm	28 cm
<b>Weight</b>	17 kg	18 kg
<b>Maximum capacity</b>	120 kg	120 k



**Aluminum Wheelchair**

Art. #VCWK9AH **PRICE: 30% below price from the market**

**PRODUCT CHARACTERISTICS:**

- foldable, light aluminium frame (double cross construction),
- nylon black upholstery,
- anti-bedsore cushion,
- lifted armrests,
- lifted and removable footrests,
- double brake system (for patient and caregiver),
- safety belt,
- pneumatic quick-release rear wheels 24",
- full front wheels 200x45,
- rear pocket,
- frame colour: platinum, blue, green.

**USE:**

General use wheelchair, at home and for strolls.

	16'	18'
<b>Seat width</b>	41 cm	45 cm
<b>Seat depth</b>	44 cm	44 cm
<b>Wheelchair height</b>	93 cm	93 cm
<b>Wheelchair width</b>	62 cm	65 cm
<b>Wheelchair length</b>	79 cm	79 cm
<b>Width of folded wheelchair</b>	32 cm	32 cm
<b>Weight</b>	16 kg	17 kg
<b>Maximum capacity</b>	100 kg	100 kg



**PRODUCT CHARACTERISTICS:**

- light, aluminium frame (double cross construction),
- soft, nylon seat and backrest,
- soft cushion,
- lifted, removable footrests of adjustable length,
- footrests with calf strap,
- lifted, removable armrests,
- safety belt,
- full, front wheels 200x50,
- pneumatic quick-release rear wheels 24",
- aluminium handrims,
- double brake system (for patient and caregiver),
- frame colour: platinum.

**USE:**

General use wheelchair, at home and for strolls.



	16"	18"	50"
<b>Seat width</b>	40 cm	45 cm	50 cm
<b>Seat depth</b>	40 cm	40 cm	40 cm
<b>Wheelchair height</b>	92 cm	92 cm	92 cm
<b>Wheelchair width</b>	55 cm	60 cm	65 cm
<b>Wheelchair length</b>	80 cm	80 cm	80 cm
<b>Width of folded wheelchair</b>	28 cm	28 cm	28 cm
<b>Weight</b>	16 kg	16 kg	16 kg
<b>Maximum capacity</b>	140 kg	140 kg	140 kg



#### PRODUCT CHARACTERISTICS:

- foldable, aluminium frame (double cross construction),
- removable black and red, or black and blue upholstery which absorbs perspiration,
- Velcro® straps allow to adjust backrest depth,
- anti-bedsore cushion,
- lifted armrests,
- lifted and removable footrests,
- calf strap,
- double brake system (for patient and caregiver),
- safety belts,
- pneumatic quick-release rear wheels 24",
- full front wheels 200x45,
- aluminium handrims,
- rear pocket,
- anti-flip wheels,
- back folding mechanism situated next to armrests allows to reduce wheelchair size during transport (VCWK9AF).

#### USE:

General use wheelchair, at home and for strolls.

	16"	18"
<b>Seat width</b>	41 cm	45 cm
<b>Seat depth</b>	42 cm	42 cm
<b>Wheelchair height</b>	90 cm	90 cm
<b>Wheelchair width</b>	61 cm	66 cm
<b>Wheelchair length</b>	81 cm	81 cm
<b>Width of folded wheelchair</b>	34 cm	34 cm
<b>Weight</b>	16 kg	17 kg
<b>Maximum capacity</b>	100 kg	100 kg



#### PRODUCT CHARACTERISTICS:

- foldable, aluminium frame of cross construction,
- detachable nylon black upholstery,
- Velcro® straps allow to adjust backrest depth,
- thick, soft cushion,
- safety belt,
- lifted, removable, angle adjustable footrests and calf strap,
- adjustable footrest length,
- adaptable position of rear wheel axle,
- lifted, detachable armrests,
- adjustable handle height,
- rear barrier which provides additional stabilization and facilitates folding the wheelchair,
- rear anti-tip wheels of adjustable height,
- full front wheels, size 200x50,
- pneumatic quick-release rear wheels 24",
- pump,
- cane holder,

- brakes of adjustable strength,
- aluminium handrims.

**USE:**

General use wheelchair, at home and for strolls. Suitable for patients after leg (legs) amputation, for patients during post-operative period.

<b>Seat width</b>	46 cm
<b>Seat depth</b>	43 cm
<b>Wheelchair width</b>	65 cm
<b>Wheelchair height</b>	103,5 cm
<b>Wheelchair length</b>	109 cm
<b>Weight</b>	19,9 kg
<b>Maximum capacity</b>	115 kg



**Aluminum Wheelchair ActiveSport**

Art. # VCWK9AS **PRICE: 30% below price from the market**

**PRODUCT CHARACTERISTICS:**

- light, foldable, aluminium frame of cross construction,
- black, detachable nylon upholstery,
- Velcro® straps allow to adjust backrest depth,
- foldable backrest,
- corrugated anti-bedsore cushion,
- safety belt,
- swing-away, detachable footrests with calf strap,
- adjustable footrest length,
- adjustable wheel position and height,
- height adjustable rear anti-tip wheels,
- full front 6" wheels,
- pneumatic quick-release rear wheels 24",
- pump,
- brakes of adjustable pressure,

- aluminium handrims,
- modern construction and design,
- frame colour: silver, blue.

**USE:**

General use comfortable wheelchair for people with reduced mobility, to use at home and for strolls.

	15"	16"	18"
Seat width	38 cm	40 cm	45 cm
<b>Seat depth</b>	39 cm	39 cm	39 cm
<b>Wheelchair height</b>	92 cm	92 cm	92 cm
<b>Wheelchair width</b>	62 cm	64 cm	70 cm
<b>Wheelchair length</b>	78 cm	78 cm	78 cm
<b>Width of folded wheelchair</b>	33 cm	33 cm	33 cm
<b>Weight</b>	14 kg	14,5 kg	15 kg
<b>Maximum capacity</b>	115 kg	115 kg	115 kg



**Tornado Active Wheelchair (Standard Version)**

Art. # MTBWA1 **PRICE: 30% below price from the market**





Aviator Active Wheelchair

Art. # MTBWB1 PRICE: 30% below price from the market

**PRODUCT CHARACTERISTICS:**

- smooth control of the gravity centre
- ergonomic seat (flat near the back, then inclined)
- narrower front frame part allows for keeping your legs closer together,

**PRODUCT CHARACTERISTICS:**

- step control of the gravity centre
- inclined seat,
- wider front frame part allows for keeping your legs wider apart,

<b>Seat width</b>	36,	38,	40,	42,	44 cm
<b>Backrest height</b>	25,	28,	31,		35 cm
<b>Seat incline angle</b>	7°,	9°,	11°,	13°	
<b>Rear wheel camber angle</b>	0°,	1°,	2°,	3°	

**COMMON FEATURES:**

- lightweight aluminium frame,
- back folds level with the seat,
- adjustable back angle (-4°, -3°, 0°, 3°, 5°, 7°),
- seat and back upholstery made of spacer fabric which facilitates air flow, fixed with Velcro® tapes to allow seat and back adjustment,
- soft padded push handles, calf protection belt and shin section of the frame,
- fibreglass side panels,

- smooth footrest adjustment,
- four footrest depths available: 10, 12, 14, 16 cm (selectable)
- inflatable rear wheels with a diameter of 24", mounted on quick-release axles to allow quick assembly and disassembly, 36 spokes, reinforced rims, grey and black Schwalbe Rightrun tyres,
- coated aluminium handrims,
- 4" front caster wheels, aluminium rim, rubber width 34 or 24 mm,
- adjustable front wheel height,
- adjustable brake position (vertical or horizontal)
- frame colour: white, black, matte black, silver, champagne (free customisation).

#### OPTIONAL ELEMENTS:

- front wheels,
- fork,
- handrims
- handrim covers,
- transparent spokes cover,
- headrest,
- anti pressure sore cushion,
- footplate,
- side guards,
- breaks,
- handles,
- safety belt,
- alcantara or leather upholstery,
- anti-tip system,
- frame colour.



**Special Wheelchair Stabilizing Back and Head with**

**Toilet Function**

Art. # VCWK7T **PRICE: 30% below price from the market**

## PRODUCT CHARACTERISTICS:

- recliner type wheelchair, with the possibility of sloping the back from sitting to horizontal position,
- foldable steel frame (double cross construction),
- black seat made of nylon or polyurethane material,
- divided seat with demountable element, which enables to get sanitary hole, removable, lifted, angle adjustable footrests with calf pads,
- removable armrests,
- removable upper part of backrest together with headrest,
- anti-tip system,
- full front wheels 200x45,
- pneumatic quick-release rear wheels 24",
- aluminium handrims,
- double brake system (for patient and caregiver),
- sanitary bucket,
- tray,
- pump.

## USE:

General use wheelchair, to use at home and for strolls with toilet function.

<b>Seat width</b>	46 cm
<b>Seat depth</b>	46 cm
<b>Wheelchair width</b>	69 cm
<b>Wheelchair height</b>	120 cm
<b>Wheelchair length</b>	136 cm
<b>Width of folded wheelchair</b>	32 cm
<b>Weight</b>	28,5 kg
<b>Maximum capacity</b>	90 kg



Special WheelChair Stabilizing Back and Head

Art. # VCWK703 PRICE: 30% below price from the market

#### PRODUCT CHARACTERISTICS:

- recliner type wheelchair, with the possibility of sloping the back from sitting to horizontal position,
- foldable aluminium frame (double cross construction),
- black seat made of nylon with thick, soft cushion,
- safe belt,
- removable, lifted, angle and length adjustable footrests with calf pads,
- removable armrests,
- removable upper part of backrest together with headrest,
- rear height adjustable anti-tip wheels,
- pneumatic quick-release rear wheels 24",
- full front wheels 200x45,
- aluminium handrims,
- double brake system (for patient and caregiver), equipped with brake pads,
- rear barrier which provides additional stabilization,
- tray,
- pump.



**USE:**General use wheelchair for people with reduced mobility, to use at home and for strolls.

<b>Seat width</b>	46 cm
<b>Seat depth</b>	46 cm
<b>Wheelchair height</b>	121 cm
<b>Wheelchair height without upper part</b>	94 cm
<b>Backrest height</b>	42cm
<b>Wheel chair length</b>	135 - 168 cm
<b>Weight ~</b>	23 kg
<b>Maximum capacity</b>	100 kg



**Special Wheelchair Stabilizing Back and Head**

Art. # VCWK702 PRICE: 30% below price from the market

## PRODUCT CHARACTERISTICS:

- recliner type wheelchair, with the possibility of sloping the back from sitting to horizontal position,
- foldable aluminium frame of cross construction,
- black seat made of nylon with thick, soft cushion,
- safe belt,
- removable, lifted, angle and length adjustable footrests with calf pads,
- removable armrests,
- removable upper part of backrest together with headrest,
- rear height adjustable anti-tip wheels,
- pneumatic rear 24" quick-release wheels,
- full front wheels 8",
- aluminium handrims,
- double brake system (for patient and caregiver), equipped with brake pads,
- rear barrier which provides additional stabilization,
- pump.

### USE:

General use wheelchair for people with reduced mobility, to use at home and for strolls.

<b>Seat width</b>	44 cm
<b>Seat depth</b>	40 cm
<b>Wheelchair height</b>	122 cm
<b>Wheelchair width</b>	68 cm
<b>Wheelchair length</b>	114 – 198 cm
<b>Width of folded wheelchair</b>	35 cm
<b>Weight</b>	27 kg
<b>Maximum capacity</b>	110 kg



**Special Wheelchair Stabilizing Back and Head**

Art. # VCWK7C PRICE: 30% below price from the market

#### **PRODUCT CHARACTERISTICS:**

- foldable aluminium frame with lateral reinforcement (cross construction),
- detachable, anti-bedsore upholstery made of memory foam equipped with waterproof cover,
- Velcro® straps allow to adjust seat and backrest depth,
- two mechanisms of angle adjustable backrest – for a patient and for a caregiver,
- adjustable wheelbase,
- safety belt,
- lifted, detachable, angle adjustable footrests and calf pads,
- adjustable position and height of rear wheels,
- height adjustable footrests,
- adjustable height of front wheels axle,
- detachable, height adjustable armrests made of soft, comfortable and easy to wash fabric,
- smooth regulation of armrest horizontal position,
- rear barrier which provides additional stabilization and facilitates folding the wheelchair,
- bottom barrier protects the wheelchair from folding during use,
- height adjustable handles,
- detachable height and angle adjustable headrest,
- brakes of adjustable strength,

- aluminium handrims,
- full front wheels 200x45,
- pneumatic quick-release rear wheels 24",
- rear, height adjustable anti-tip wheels,
- pump,
- tray.

**USE:**

General use comfortable and versatile wheelchair for people with reduced mobility, to use at home and for strolls.

	16"	18"
<b>Seat width</b>	42 cm	46 cm
<b>Seat depth</b>	50 cm	50 cm
<b>Wheelchair height</b>	100,8 cm	100,8 cm
<b>Wheelchair width</b>	65 cm	71,5 cm
<b>Wheelchair length</b>	114 cm	114 cm
<b>Width of folded wheelchair</b>	36 cm	36 cm
<b>Weight</b>	28.7 kg	28.7 kg
<b>Maximum capacity</b>	136 kg	136 kg





## PRODUCT CHARACTERISTICS:

- stable, reinforced steel frame,
- foldable backrest,
- detachable, anti-bedsore upholstery made of 'memory foam' equipped with waterproof cover,
- anatomically moulded backrest and seat,
- smooth adjustment of backrest height,
- smooth adjustment of seat angle (together with backrest),
- smooth adjustment of backrest angle,
- height adjustable footrest,
- swing-away, detachable, angle adjustable footrests and calf pads,
- height adjustable front wheels axle,
- detachable, height adjustable armrests made of soft, comfortable and easy to wash fabric,
- smooth adjustment of armrest horizontal position,
- adjustable seat width,
- double handle of adjustable height,
- lateral wedges of adjustable horizontal and vertical position,
- rear, height adjustable anti-tip wheels,
- detachable, height and tilt adjustable headrest,
- brakes of adjustable strength,
- aluminium handrims,
- full front wheels, size 200x50,
- pneumatic, quick-release rear wheels 24",
- pump.

## USE:

General use comfortable and versatile wheelchair for people with reduced mobility, to use at home and for strolls.

	43	50
<b>Seat width</b>	36 - 43 cm	45 - 50 cm
<b>Seat depth</b>	42 - 48 cm	42 - 48 cm
<b>Wheelchair height</b>	96 cm	96 cm
<b>Wheelchair width</b>	54,5 cm	64,5 cm
<b>Wheelchair length</b>	108 cm	108 cm
<b>Weight</b>	42 kg	42 kg
<b>Maximum capacity</b>	136 kg	136 kg



Wheelchair

Baczus Relax Aluminium Children's

Art. # VCWR211R **PRICE: 30% below price from the market**

#### PRODUCT CHARACTERISTICS:

- foldable, aluminium frame of cross construction,
- detachable, easy to wash nylon upholstery,
- detachable seat and backrest,
- reclining backrest (together with seat),
- height and tilt adjustable headrest,
- detachable side supports,
- abduction wedge of adjustable depth,
- footrests of adjustable height, calf pads,
- side stabilizing wedges,
- swing-away, angle adjustable armrests,
- security suspenders,
- height adjustable front wheel axle,
- brakes of adjustable strength,
- aluminium handrims,
- full front 6" wheels,
- rear height adjustable anti-tip wheels,
- pneumatic quick-release rear wheels 20",
- pump,
- tray,
- frame colour: red, blue.

**USE:**For children with cerebral palsy, spina bifida, muscular diseases.

	13"	15"
<b>Seat width</b>	34 cm	38 cm
<b>Seat depth</b>	36 cm	36 cm
<b>Wheelchair width</b>	97,5 cm	97,5 cm
<b>Wheelchair length</b>	54,5 cm	64,6 cm
<b>Wheelchair height</b>	87,5 cm	87,5 cm
<b>Weight</b>	22,8 kg	22,8 kg
<b>Maximum capacity</b>	90 kg	90 kg



**Baczus Applause Children Rehabilitative**

**Stroller**

Art # VCWR211 **PRICE: 30% below price from the market**

**PRODUCT CHARACTERISTICS:**

- aluminium and steel frame,
- height and width adjustable backrest,
- angle adjustable backrest (sitting-horizontal position),
- depth and width adjustable seat,
- length and angle adjustable footrest,
- height and angle adjustable handle for caregiver,
- swivel front wheels with lock,
- front guard and positing vest,
- guard-rail and stabilising braces,

- detachable abduction wedge and headrest,
- spring brace,
- standard equipment: canopy, additional equipment: umbrella, rain cover, sleeping bag,
- available colours: grey/red.

**USE:**

For children with cerebral palsy, spina bifida, muscular diseases.

<b>Footrest height</b>	17 - 35 cm
<b>Seat width</b>	26 - 44 cm
<b>Seat depth</b>	30 - 42 cm
<b>Backrest height</b>	43 - 56 cm
<b>Range of backrest angle</b>	0° - 90°
<b>Stroller max width</b>	66 cm
<b>Stroller max length</b>	100 cm
<b>Stroller max height</b>	99 cm
<b>Dimensions when folded (length x width x height)</b>	104x66x50 cm
<b>Front wheel diameter</b>	20 cm
<b>Rear wheel diameter</b>	32 cm
<b>Stroller weight</b>	21,4 kg
<b>Maximum capacity</b>	35 kg





#### PRODUCT CHARACTERISTICS:

- aluminium, foldable frame,
- polyester fabric seat,
- polyester backrest of adjustable inclination (two positions),
- length adjustable, aluminium footrests with calf strap,
- swivel front wheels (foam or pneumatic) with possibility of direction lock,
- removable rear wheels (foam or pneumatic) with shock absorber and independent brake,
- front guard and seat belt,
- abduction wedge,
- headrest,
- standard equipment: canopy, sleeping bag, rain cover, shopping basket, bag,
- colors available: black/blue, black/green, black/orange, black/pink.

#### USE:

Children's strolling cart for children with Cerebral Palsy, Spina Bifida and muscle disorders.

<b>Footrest height</b>	20 - 35 cm
<b>Seat width</b>	39 cm
<b>Seat depth</b>	30 cm
<b>Backrest height</b>	54 cm
<b>Range of backrest angle</b>	0 - 40°
<b>Max stroller width</b>	60 cm
<b>Max stroller length</b>	100 cm
<b>Max stroller height</b>	116 cm
<b>Dimension while folded</b>	126x50x40 cm
<b>Front wheels diameter</b>	18.5 cm
<b>Rear wheels diameter</b>	22.2 cm
<b>Wheelchair weight</b>	11.5 kg
<b>Maximum capacity</b>	30 kg



**Umbrella 2 Children Rehabilitative Stroller**

**PRICE: 30% below price from the market**



**Junior Children Rehabilitative Stroller**

Art. # DRVGOE **PRICE: 30% below price from the market**

#### **PRODUCT CHARACTERISTICS:**

- foldable, steel and aluminium frame,
- height and width adjustable backrest,
- angle adjustable backrest (sitting-horizontal position),
- adjustable seat width and depth,
- width adjustable lateral supports,
- reversible seat unit,
- length and position adjustable footrest,
- adjustable handle,
- front foam swivel wheels with quick-release system,
- rear foam wheels with shock absorbers, lock and quick-release system,
- height adjustable headrest,
- front guard and positioning vest,
- abduction wedge,
- standard equipment: canopy, sleeping bag, rain cover, shopping basket,
- available colours: black/gray, black/green/ black/blue.

#### **USE:**

Children's strolling cart for children with cerebral palsy, spinablifidawith myelocole, muscle dysfunctions, lower limb dysfuctions.

Footrest height	31 - 40 cm
Seat width	25 - 30 cm
Seat depth	27 - 33 cm
Backrest height	47 - 51 cm
Range of backrest angle	0° - 90°
Max stroller width	66 cm
Max stroller length	120 cm
Max stroller height	117 cm
Dimensions when folded (length x width x height)	110x66x65 cm
Front wheel diameter	22 cm
Rear wheel diameter	27 cm
Stroller weight	18,5 kg
Maximum capacity	30 kg





## PRODUCT CHARACTERISTICS:

- wooden rehabilitative chair,
- angle and depth adjustable seat,
- adjustable lateral supports,
- angle adjustable backrest (sitting-horizontal position),
- angle and length adjustable footrests,
- angle adjustable foot platform with heel support and safety belt,
- adjustable abduction wedge,
- height and angle adjustable headrest,
- height and angle adjustable tray,
- positioning vest,
- swivel wheels with lock,
- adjustable headrest,
- height adjustable support of the cervical spine,
- shelf under the seat.

## USE:

Rehabilitative chair for disabled children, children diagnosed with nervous system disorder, for children who have difficulty sitting.

<b>Seat height (with headrest)</b>	46 - 62 cm
<b>Seat depth</b>	25 - 31 cm
<b>Footrest height</b>	15 - 31 cm
<b>Armrest height</b>	68 - 76 cm
<b>Seat width</b>	15 - 35 cm
<b>Backrest width</b>	15 - 27 cm
<b>Chair width</b>	58 cm
<b>Chair length</b>	95 - 150 cm
<b>Chair weight</b>	27 kg
<b>Maximum capacity</b>	35 kg